

AllisonAndAssociates.com  
 Ph: 501-753-4648  
 Toll Free: 800-264-0648  
 Fax: 501-753-6062  
 Cell: 501-690-8660



# Allison and Associates Manufacturing Company Inc

## APPROVED CLEANERS

**Important Notes on Last Page!**

MANUFACTURER	PRODUCT	INSTRUCTIONS FOR USE WITH CARRFLEX® EXTREME	INSTRUCTIONS FOR USE WITH STANDARD CARRFLEX®
Ecolab®	Actichlor Plus	0.2% Actichlor Plus (0.2:100 Water)	0.2% Actichlor Plus (2000PPM)
Quantum Pacific Limited	Actimax	According to Manufacturer Spec	According to Manufacturer Spec
ABC Compounding Co., Inc.	Aero Stain Away	According to Manufacturer Spec	*
Milton®	Anti-Bacterial Tablets	1 Tablet:2L Water	1 Tablet:2L Water
Ecolab®	Asepti Active Liquid	According to Manufacturer Spec	According to Manufacturer Spec
Ecolab®	Asepti-Wipe II	Ready to Use	Ready to Use
Synergy Health	Azomax™ Multi Surface Wipes	According to Manufacturer Spec	Ready to Use
Synergy Health	Azowipe™ Hard Surface Wipes	According to Manufacturer Spec	Ready to Use
Clorox®	Bleach Wipes (Sodium Hypochlorite 0.1 - 1%)	According to Manufacturer Spec	*
Spartan Chemical Company Inc.	BNC-15™	According to Manufacturer Spec	According to Manufacturer Spec
Ameriderm Laboratories Limited	Care Source Disinfectant	According to Manufacturer Spec	*
Metrex® Research	CaviCide®	According to Manufacturer Spec	According to Manufacturer Spec
Metrex® Research	CaviCide1™	According to Manufacturer Spec	According to Manufacturer Spec
Helix Solutions	Chlor-Clean	1000PPM (1 Tablet per 1 Litre Water)	*
Agar Cleaning Systems Pty. Ltd.	Chloradet	0.6% Chloradet (1:160 Water)	0.1% Chloradet (1:1000 Water)
Ecolab®	Chlorwhite	1% Chlorwhite (1:100 Water)	*
Appollo Corporation	CID A-L™ II Quaternary Disinfectant	According to Manufacturer Spec	According to Manufacturer Spec
Clorox Disinfecting Wipes	Clorox Disinfecting Wipes	Ready to Use	According to Manufacturer Spec
Alco Mandala Products	DC-7	According to Manufacturer Spec	According to Manufacturer Spec
Delaval	Delaval Acid Wash	According to Manufacturer Spec	*
Clinell®	Detergent Wipes	According to Manufacturer Spec	According to Manufacturer Spec
ArjoHuntleigh	Disinfectant Cleanser IV	According to Manufacturer Spec	According to Manufacturer Spec
Clorox®	Dispatch® Disinfectant	According to Manufacturer Spec	According to Manufacturer Spec
Tristel Solutions LTD	Distel Surface Disinfectant (Formerly Trigene®)	According to Manufacturer Spec	According to Manufacturer Spec
Sealed Air Diversey Care	Divercleanse	3.1% Divercleanse (3.1:100 Water)	3.1% Divercleanse (3.1:100 Water)
Buckeye® International, Inc.	Eco E23 S23	256 to 1 dilution with water	256 to 1 dilution with water
Velocity Chemicals	Eco-Pine	5% Ecopine (1:20 Water)	5% Ecopine (1:20 Water)
Clorox HEALTHCARE®	Fuzion™	According to Manufacturer Spec	According to Manufacturer Spec
Spartan Chemical Company Inc.	GS Neutral Disinfectant #3502	According to Manufacturer Spec	According to Manufacturer Spec
Clinell®	Hand & Surface Sanitizer	According to Manufacturer Spec	According to Manufacturer Spec
Spartan Chemical Company Inc.	hdq®C2	According to Manufacturer Spec	According to Manufacturer Spec
Brighton Professional	Hepastat 256	According to Manufacturer Spec	*

AllisonAndAssociates.com  
 Ph: 501-753-4648  
 Toll Free: 800-264-0648  
 Fax: 501-753-6062  
 Cell: 501-690-8660



# Allison and Associates

## Manufacturing Company Inc

### APPROVED CLEANERS

**Important Notes on Last Page!**

MANUFACTURER	PRODUCT	INSTRUCTIONS FOR USE WITH CARRFLEX® EXTREME	INSTRUCTIONS FOR USE WITH STANDARD CARRFLEX®
(Various)	Hydrogen Peroxide 5%	Ready to Use	*
Appollo Corporation	Hygena™	According to Manufacturer Spec	According to Manufacturer Spec
Ecolab®	Incidin® Extra	According to Manufacturer Spec	*
Ecolab®	Incidin® Foam	According to Manufacturer Spec	*
Ecolab®	Incidin® Plus	According to Manufacturer Spec	*
(Various)	Isopropyl Alcohol 70%	Ready to Use	*
Kimberly-Clark®	Isowipes®	Ready to Use	Ready to Use
Schulke & Mayr GMBH	Kodan® Tinkture Forte	According to Manufacturer Spec	*
Legge Systems	Legcide HB Plus	According to Manufacturer Spec	According to Manufacturer Spec
Whiteley Medical	Lemex	According to Manufacturer Spec	According to Manufacturer Spec
Livingstone International (PTY) LTD	LIV-Wipes Jumbo	According to Manufacturer Spec	Ready to Use
Sulco LTD	Mediwipes	According to Manufacturer Spec	Ready to Use
Pro Restore® Products	Microban X-590 Institutional Spray	According to Manufacturer Spec	According to Manufacturer Spec
Medline Industries Inc.	Micro-Kill™ AF2 Disinfecting, Deodorizing, Cleaning Wipes	Ready to Use	Ready to Use
Medline Industries Inc.	Micro-Kill™ Q3 Concentrated Disinfectant, Cleaner & Deodorizer	According to Manufacturer Spec	According to Manufacturer Spec
Anti-Slip Anywhere	Multi Kleen	According to Manufacturer Spec	According to Manufacturer Spec
Ecolab® Cleantec	Nature Clean	1: 38 Water	*
Certo®	Non Acid Bathroom Cleaner	According to Manufacturer Spec	According to Manufacturer Spec
Ecolab®	Oasis Pro 12 Neutral Cleaner	5% Oasis Pro 12 (1:20 Water)	5% Oasis Pro 12 (1:20 Water)
Schülke & Mayr GmbH	Ocetenisept®	According to Manufacturer Spec	*
Whiteley Medical	Odacon Spray	According to Manufacturer Spec	According to Manufacturer Spec
OdorXit Naturals	Odor Eliminator Spray	Ready to Use	Ready to Use
Sealed Air Diversey Care	Oxivir® TB Liquid and Wipes	According to Manufacturer Spec	*
Ecolab®	Oxycide	According to Manufacturer Spec	*
Pal® International	PAL® Detergent Wipes	According to Manufacturer Spec	Ready to Use
Pal® International	PAL® Disinfectant Wipes	According to Manufacturer Spec	Ready to Use
PCS Process Cleaning Solutions	PCS 7000	According to Manufacturer Spec	According to Manufacturer Spec
PCS Process Cleaning Solutions	PCS Spray & Wipe Hand and Surface Cleaner	Ready to Use	Ready to Use
Pure Green, LLC	Puregreen 24 (Pure Green)	According to Manufacturer Spec	According to Manufacturer Spec
Reynard Health Supplies	Reynard Neutral Detergent Wipes	According to Manufacturer Spec	According to Manufacturer Spec
Reynard Health Supplies	Reynard Surface Sanitizing Wipes	According to Manufacturer Spec	According to Manufacturer Spec
Sealed Air Diversey Care	Room Care R2L	According to Manufacturer Spec	According to Manufacturer Spec

AllisonAndAssociates.com  
 Ph: 501-753-4648  
 Toll Free: 800-264-0648  
 Fax: 501-753-6062  
 Cell: 501-690-8660



# Allison and Associates Manufacturing Company Inc

## APPROVED CLEANERS Important Notes on Last Page!

MANUFACTURER	PRODUCT	INSTRUCTIONS FOR USE WITH CARRFLEX® EXTREME	INSTRUCTIONS FOR USE WITH STANDARD CARRFLEX®
STERI-7®	S-7 Xtra Wipes	Ready to Use	Ready to Use
Buckeye® International, Inc.	Sanicare TBX™	According to Manufacturer Spec	*
Buckeye® International, Inc.	Sanicare™ Disinfectant Wipes	According to Manufacturer Spec	*
Alco Mandala Products	Sanni-Rinse	According to Manufacturer Spec	According to Manufacturer Spec
(Various)	Sodium Hypochlorite	1.3% (1.3:100 Water) (Approved Colors)	*
Hu-Friedy Mfg. Co., LLC.	Sporeclear RTU Spray or Wipes	Diluted at 8% Into Water	*
Clinell®	Sporicidal Wipes	According to Manufacturer Spec	According to Manufacturer Spec
Ecolab®	Spring Clean	5% Spring Clean (1:20 Water)	5% Spring Clean (1:20 Water)
Noble Pine Products Co.	Steri-Fab®	According to Manufacturer Spec	According to Manufacturer Spec
Hartmann	Sterillium®	According to Manufacturer Spec	*
sBioMed, LLC	Steriplex SD Disinfectant Sporicide Wipes (Without Bleach)	According to Manufacturer Spec	According to Manufacturer Spec
Sterri-Matt®	Sterri-Matt® Foam	Ready to Use	Ready to Use
Sealed Air Diversey Care	Stride® Citrus Smartdose	1:750 Water	*
Sealed Air Diversey Care	Stride® HC (JF)	According to Manufacturer Spec	*
Sealed Air Diversey Care	Suma San Conc D4A	According to Manufacturer Spec	*
Sealed Air Diversey Care	Suma San Conc D4A	According to Manufacturer Spec	*
Spartan Chemical Company Inc.	Super HDQ Neutral®	256 to 1 dilution with water	256 to 1 dilution with water
PDI	Super Sani Cloth®	According to Manufacturer Spec	Ready to Use
Anios Laboratories	Surfa 'Safe	According to Manufacturer Spec	*
Bioscience Fresche®	Surface Plus	According to Manufacturer Spec	*
Reynard Health Supplies	Surface Sanitising Wipes	Ready to Use	Ready to Use
3M	TB Quat Disinfectant	Ready to Use	Ready to Use
Buckeye® International, Inc.	Terminator™ One Step Disinfectant	According to Manufacturer Spec	*
Appollo Corporation	Therasol™	According to Manufacturer Spec	According to Manufacturer Spec
Hunters Products	TRI-DET	1:80 Water	*
(Various)	Trix (Soap)	Ready to Use	Ready to Use
Verna Care	Tuffie 5 Cleaning Wipes	According to Manufacturer Spec	Ready to Use
Apollo Corporation	Turbo Clean™	According to Manufacturer Spec	According to Manufacturer Spec
Clinell®	Universal Disinfectant Wipes	According to Manufacturer Spec	According to Manufacturer Spec
True Blue Chemicals	Use All	According to Manufacturer Spec	According to Manufacturer Spec
Whiteley Medical	Viraclean®	According to Manufacturer Spec	According to Manufacturer Spec
Johnson Diversey	Virex II 256 (US)	According to Manufacturer Spec	According to Manufacturer Spec

AllisonAndAssociates.com  
 Ph: 501-753-4648  
 Toll Free: 800-264-0648  
 Fax: 501-753-6062  
 Cell: 501-690-8660



# Allison and Associates

## Manufacturing Company Inc

### APPROVED CLEANERS

**Important Notes on Last Page!**

MANUFACTURER	PRODUCT	INSTRUCTIONS FOR USE WITH CARRFLEX® EXTREME	INSTRUCTIONS FOR USE WITH STANDARD CARRFLEX®
Dupont™	Virkon®	5% Virkon (1:20 Water)	1% Virkon (1:100 Water)
Wexford Labs	Wexcide 128	Minimum 128 to 1 dilution with water	Minimum 128 to 1 dilution with water

\*The cleaner listed is recommended for use only with CARRFLEX® EXTREME

#### IMPORTANT INFORMATION:

Information contained here is purely related to the compatibility of cleaners with CARRFLEX® materials and in no way relates to the efficacy of any cleaner for infection control or the safety of any cleaner for skin contact.

Testing info is specific to the products or formulations listed and in no way constitutes approval for other products from any brand listed.

All approvals are based on manufacturer recommended dilution with water unless otherwise specified and approval is void for any improper dilution.

If a particular cleaner/disinfectant is not listed, we advise to test the specific product in advance or to contact us to conduct testing.

#### TIPS FOR EXTENDING THE LIFE OF ANY MATTRESS SURFACE:

Never use an abrasive cloth, sponge or other similar device on the surface.

Wipe off cleaning product residue with a clean, soft cloth and rinse with fresh warm water immediately after use.

Allow the surface to dry completely before use.